PSD 1

pH MÈTRE / pH METER

Measure up

Français

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English

Thank you for purchasing your **PSD 1 pH Meter**.

For best results with your device:

- **Read** these operating instructions carefully,
- **Observe** the precautions of use.

MEANINGS OF THE SYMBOLS USED

	WARNING, risk of DANGER! The operator must refer to this user's manual whenever this danger symbol appears.
	The CE marking indicates conformity with European directives, in particular LVD and EMC.
	The rubbish bin with a line through it indicates that, in the European Union, the product must undergo selective disposal in compliance with Directive WEEE 2002/96/EC. This equipment must not be treated as household waste.

CONTENTS

1. TECHNICAL CHARACTERISTICS.....	14
2. DESCRIPTION	15
2.1 TOP.....	15
2.2 REAR PANEL	16
2.3 BOTTOM	16
3. USE.....	17
3.1 CONNECTION AND STARTING UP	17
3.2 ADJUSTMENTS AND MEASUREMENTS OF THE pH.....	17
3.2.1 <i>MANUAL TEMPERATURE CORRECTION</i>	17
3.2.2 <i>CALIBRATION AND ADJUSTMENT OF THE pH SCALE</i>	17
3.3 RECORDER OUTPUT (0-5V).....	18
4. VARIOUS POSSIBLE PROBLEMS	18
5. SAFETY RULES	19
6. MAINTENANCE.....	20
6.1 CLEANING	20
6.2 METROLOGICAL CHECK.....	20
6.3 REPAIRS	20
7. TO ORDER.....	21
7.1 DELIVERY CONDITION.....	21

1. TECHNICAL CHARACTERISTICS

2000-point liquid crystal display	
Range	0.00 to 14.00 pH
Accuracy	± 1%
Manual temperature compensation	0 to 100°C
Temperature of use	0-50°C
Maximum altitude	2000m
Humidity	5 to 80% at 35°C
Pollution class	Category 2
Category of use	2
COMPLIANCE WITH STANDARDS:	
Safety	With the power supply provided: EN 61010-1. Normal environmental conditions (indoor use)
CEM	EN 61326, class B instrument, intermittent operation.
Consumption	0.5W
Dimensions	155 x 90 x 41 mm
Weight	200g


The PSD 1 is provided with a power supply unit:

Power supply	230V 50-60Hz
Output	9V 


The PSD 1 can operate only on DC.

2. DESCRIPTION

2.1 TOP


2.2 REAR PANEL

2.3 BOTTOM

- Identification label

3. USE

3.1 CONNECTION AND STARTING UP

Unpack the instrument, the power supply unit and the operating manual.
Connect the pH electrode to the BNC connector on the back of the instrument.
Connect the power supply unit to the "power supply" connector.
(The PSD1 must be used only with the power supply unit provided).

3.2 ADJUSTMENTS AND MEASUREMENTS OF THE pH

3.2.1 MANUAL TEMPERATURE CORRECTION

Note the temperature of the sample to be analyzed.

Set the mark on the "Temperature" knob to the temperature of the medium to be analyzed or of the buffer solution.

Remark: Ideally, the buffer solutions used to calibrate the measurement system (pH meter-electrode) should be at this same temperature.

3.2.2 CALIBRATION AND ADJUSTMENT OF THE pH SCALE

3.2.2.1 Checks/precautions concerning the electrodes:

Before each pH measurement operation, make sure that:

- the glass bulb of the electrode does not contain any air bubble,
- the liquid electrolyte reference electrode is correctly filled with a potassium chloride solution of the right concentration (KCl , 1 mol.L^{-1}),
- the plug of the filling orifice of the reference electrode is open,
- the electrode is correctly connected.

3.2.2.2 Calibration:

Example: Use of pH4 and pH7 buffer solutions.

- Rinse the electrode in distilled water.
- Immerse the electrode in the pH7 buffer solution (always first).
- When the signal stabilizes, use the "Offset" knob to adjust the value displayed to 7.00
- Rinse the electrode in distilled water.
- Immerse the electrode in a pH4 buffer solution.

English

- When the signal stabilizes, use the "Slope" knob to adjust the reading to 4.00
- Rinse the electrode in distilled water.
- Immerse the electrode in the medium to be analyzed.

3.3 RECORDER OUTPUT (0-5V)

The recorder output is on the back of the instrument.

Black (-) and red (+) terminals for insulated 4mm banana plugs.

Correspondence:

pH	0-14 pH
Recorder output	0-5000mV (0-5V)

4. VARIOUS POSSIBLE PROBLEMS

PROBLEM	ADVICE
No display	Check the 230V supply. Check the connection of the power supply unit to the power connector.
The pH fails to stabilize in all solutions measured.	Check that the pH electrode is correctly immersed (reference junction, glass bulb). Check the connection of this electrode.
The instrument indicates 7 whatever the pH.	Check the condition of the pH electrode (glass bulb). Replace the pH electrode.
During calibration in the pH7 buffer solution, impossible to adjust the display to 7.00 Offset potentiometer at stop.	Check the level of the reference electrolyte. Replace the electrolyte of the 1 mol.L ⁻¹ KCl reference electrode (if liquid electrolyte).
During calibration in the pH4 or pH9 buffer solution, impossible to adjust the display to 4.00 and 9.00, respectively. Slope potentiometer at stop.	Check the position of the temperature potentiometer. Clean the electrode by dipping it a few times in a dilute hydrochloric acid solution (0.1 mol.L ⁻¹ HCl).

5. SAFETY RULES

- The mains supply must have the following characteristics: 230V \pm 10% 50-60Hz-5W.
- The interior of the instrument must always be kept clean and dry.
- If the instrument is used in a way not in conformity with the specifications, the protection provided by the instrument may be impaired.
- Disconnect the instrument before opening the housing.
- The power supply unit serves as voltage disconnect device.
- The PSD 1 pH-meter is designed for people familiar with good laboratory practice. If the PSD 1 is not used in accordance with these operating instructions, the protection provided by the equipment may be impaired.
- The operation of the instrument can be perturbed by radiated electric fields or conducted emissions.
- The instrument must not be used in a context requiring permanent operation without human supervision.
- This instrument produces, uses, and can emit RF energy, and if not installed and used in accordance with the operating instructions, it can interfere with radio communications.
- In a residential area, the user of this equipment will probably cause interference, in which case the user must, at his/her own expense, do everything necessary to remedy the interference.

6. MAINTENANCE

The instrument contains no parts that can be replaced other than by trained and accredited personnel. Any unauthorized repair or replacement of a part by an "equivalent" may gravely impair safety.

6.1 CLEANING

Use a cloth moistened with clean water or a neutral detergent to wipe the instrument, then wipe with a dry cloth.

Do not use the instrument again until it is completely dry.

6.2 METROLOGICAL CHECK

Like all measuring or testing devices, the instrument must be checked regularly.

This instrument should be checked at least once a year. For checking and calibration, contact one of our accredited metrology laboratories (information and contact details available on request), at our Chauvin Arnoux subsidiary or the branch in your country.

6.3 REPAIRS

For all repairs before or after expiry of warranty, please return the device to your distributor.

7. TO ORDER

7.1 DELIVERY CONDITION

- 1 PSD 1 pH Meter
- 1 user manual
- 2 pH buffer solutions
- 1 9V power supply unit

All packed in cardboard box.

FRANCE

Chauvin Arnoux Group

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75876 PARIS Cedex 18

Tél : +33 1 44 85 44 85

Fax : +33 1 46 27 73 89

info@chauvin-arnoux.com

www.chauvin-arnoux.com

INTERNATIONAL

Chauvin Arnoux Group

Tél : +33 1 44 85 44 38

Fax : +33 1 46 27 95 69

Our international contacts

www.chauvin-arnoux.com/contacts