VASD 40

POTENTIOSTAT

Measure up

Français

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Thank you for purchasing your **Potensiostat VASD 40**.

For best results with your device:

- **Read** these operating instructions carefully,
- **Observe** the precautions of use.

MEANINGS OF THE SYMBOLS USED

	<p>WARNING, risk of DANGER! The operator must refer to this user's manual whenever this danger symbol appears.</p>
	<p>The CE marking indicates conformity with European directives, in particular LVD and EMC.</p>
	<p>The rubbish bin with a line through it indicates that, in the European Union, the product must undergo selective disposal in compliance with Directive WEEE 2002/96/EC. This equipment must not be treated as household waste.</p>

CONTENTS

1. GENERAL INSTRUCTIONS	17
1.1 SAFETY INSTRUCTIONS	17
1.1.1 <i>BEFORE USE</i>	17
1.1.2 <i>DURING USE</i>	17
2. TECHNICAL CHARACTERISTICS.....	18
3. PRINCIPLE OF OPERATION.....	19
4. DESCRIPTION OF THE INSTRUMENT.....	20
4.1 FRONT PANEL.....	20
4.2 REAR PANEL	21
5. USE.....	22
5.1 CONNECTION AND STARTING UP	22
5.2 ADJUSTMENT	22
5.3 MEASUREMENT	22
5.3.1 <i>RECORDER OUTPUTS</i>	23
6. TESTS.....	24
7. VARIOUS POSSIBLE PROBLEMS	25
8. MAINTENANCE.....	26
8.1 CLEANING	26
8.2 METROLOGICAL CHECK.....	26
8.3 REPAIRS.....	26
9. TO ORDER.....	27
9.1 DELIVERY CONDITION.....	27

1. GENERAL INSTRUCTIONS

The potentiostat is a class B instrument per standard EN 61326-1.

(Intermittent operation)

The power supply must have the following characteristics: 230V \pm 10%, 50-60Hz.

The power supply unit serves as disconnecting switch.

Conducted interference in the power grid or powerful electromagnetic interference may cause fluctuations in the readings.

The potentiostat is intended for persons familiar with good laboratory practice.

If the instrument is not used in accordance with these operating instructions, the protection provided by the equipment may be impaired.

This instrument produces, uses, and can emit RF energy, and if it is not installed and used in accordance with the user manual, it can interfere with radio communications, in which case the user must, at its own expense, take all necessary measures to eliminate the interference.

1.1 SAFETY INSTRUCTIONS

1.1.1 BEFORE USE

The use of the VASD 40 potentiostat requires the user to observe the usual safety rules:

- To protect the use against electrical dangers,
- And to protect the instrument from operator error.

1.1.2 DURING USE

- The instrument must not be exposed to a temperature above 50°C.

2. TECHNICAL CHARACTERISTICS

2000-point liquid crystal display	
Potentials adjustable from:	-1999 to +1999 mV
Current measurement ranges	2000 μ A, 20mA
Accuracy	\pm 1%
Temperature of use	0 to 50°C
Maximum altitude	2000m
Humidity	5 to 80% at 35°C
Pollution class	Category 2
Category of use	2
Consumption	0.3W
Dimensions	204 x 150 x 53 mm
Weight	460g

The Potentiostat VASD 40 is provided with a power supply unit.

Power supply	230V – 50/60Hz
Output	9V 

The VASD 40 works in DC only.


3. PRINCIPLE OF OPERATION

The potentiostat has been developed for three-electrode amperometry. The instrument is connected to a work electrode, **T**, where there occurs a reaction with the electroactive species in the solution being investigated, to a reference electrode, **R**, and to an auxiliary electrode (or counter-electrode), **A**.

The instrument displays the potential difference \underline{E} between the work electrode and the reference electrode, or $\underline{E}(\text{work}) - E(\text{reference})$

The auxiliary electrode ensures that no current flows in the reference electrode. Thus, a voltage \underline{U} applied between the work electrode and the auxiliary electrode generates a current \underline{I} of this current is displayed for a potential $\underline{E}(\text{work})$ of the work electrode, in a given solution and according to the electrode materials used.


The user can vary the imposed voltage \underline{U} between the work electrode and the auxiliary electrode in order to observe the current \underline{I} generated over a range of potentials, positive or negative with respect to the reference.

Amperometric mounting Principle

4. DESCRIPTION OF THE INSTRUMENT


4.1 FRONT PANEL

4.2 REAR PANEL

3 insulated 4mm terminals
for electrodes

2 insulated 4mm terminals
for measurement of the potential

2 insulated 4mm terminals
for measurement of the current

Power input
(2.1mm jack)

5. USE

5.1 CONNECTION AND STARTING UP

- Unpack the instrument, the mains adapter, and the user manual.
- Connect the three electrodes to their respective connectors.
- Connect the 2.1mm plug of the mains adapter to the potentiostat.

5.2 ADJUSTMENT

The range of potentials in which the measurement will be made can be chosen to suit the manipulation: from $0 \pm 500\text{mV}$ to $0 \pm 2000\text{mV}$.

Before the range can be chosen, it is first necessary to turn the potential adjustment knob to the maximum, then use a screwdriver to adjust the screw of the potential range adjustment potentiometer so that the reading on the mV display unit is the desired maximum value.

The μA -mA range switch can be used to change ranges while making measurements, if that proves necessary.

5.3 MEASUREMENT

Plot of the anodic circuit:

- switch set to mV+
- Set $I = 0$ using the knob for adjustment of the potential difference between the work electrode and the reference electrode.
- Vary I in the direction $I > 0$ using the same knob; note the (I, U) couples.

Plot of the cathodic circuit:

- Switch set to mV+
- As before, set $I = 0$
- Vary I in the direction $I < 0$ using the same knob, U then decreases down to 0 ; to continue the measurements when $E(\text{Work}) < E(\text{Reference})$, change the switch over from mV+ to mV-

Remark: if dioxygen is evolved during the plotting of the anodic circuit (or dihydrogen, in the case of the cathodic circuit), remember to eliminate it before starting another series of measurements.

5.3.1 RECORDER OUTPUTS

The potentiostat has two 0-5V recorder outputs.



Black \ominus and red \oplus terminals for insulated 4mm banana plugs.

Correspondence:

Display		Recorder output
	Potential	
-2000mV		0V
000mV		+2.5V
+2000mV		+5V
	Current	
-2000 μ A		0V
000 μ A		+2.5V
+2000 μ A		+5V

6. TESTS

EXECUTED BY THE LYCÉE MARCELIN BERTHELOT, 94 ST MAUR DES FOSSÉS

7. VARIOUS POSSIBLE PROBLEMS

PROBLEM	ADVICE
No display	Check the 230V supply. Check the connection of the power supply unit to the power connector
The current display unit indicates 1 (value greater than 2000)	Change to the mA range.
The potential display unit indicates 1	Turn the potential adjustment knob anticlockwise.

8. MAINTENANCE

The instrument contains no parts that can be replaced other than by trained and accredited personnel. Any unauthorized repair or replacement of a part by an "equivalent" may gravely impair safety.

8.1 CLEANING

Use a cloth moistened with clean water or a neutral detergent to wipe the instrument, then wipe with a dry cloth.

Do not use the instrument again until it is completely dry.

8.2 METROLOGICAL CHECK

Like all measuring or testing devices, the instrument must be checked regularly.

This instrument should be checked at least once a year. For checking and calibration, contact one of our accredited metrology laboratories (information and contact details available on request), at our Chauvin Arnoux subsidiary or the branch in your country.

8.3 REPAIRS

For all repairs before or after expiry of warranty, please return the device to your distributor.

9. TO ORDER

9.1 DELIVERY CONDITION

- 1 VASD 40 potentiostat
- 1 user manual
- 1 9V power supply unit

All packed in a cardboard box.

FRANCE

Chauvin Arnoux Group

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75876 PARIS Cedex 18

Tél : +33 1 44 85 44 85

Fax : +33 1 46 27 73 89

info@chauvin-arnoux.com

www.chauvin-arnoux.com

INTERNATIONAL

Chauvin Arnoux Group

Tél : +33 1 44 85 44 38

Fax : +33 1 46 27 95 69

Our international contacts

www.chauvin-arnoux.com/contacts