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
Thank you for purchasing a **separate pH electrode**.

For the best results with your electrode:

- **read** this user manual carefully,
- **comply** with the precautions for use.

This electrode can be used in association with a reference electrode to measure the pH of a solution.

Certain electrodes may contain mercury.

These products are highly poisonous, even at low doses. They may have very varied physical effects: nausea, vomiting, headache, loss of consciousness or other more significant problems leading to death.

These products may be included in one or more of the following categories: carcinogenic, mutagenic, toxic for human reproduction. They may also modify the functioning of certain organs (liver, nervous system), attack the lungs and cause allergies (asthma).

These products may have harmful effects on the environment, in particular on organisms in aquatic environments: fish, crustaceans, algae and other aquatic plants.

The CE marking indicates compliance with the European directives, particularly LV and EMC.

The crossed-out dustbin means that, in the European Union, the product is the subject of selective disposal as per directive DEEE 2002/96/EC: this equipment must not be treated as household waste.

USE

- Remove the storage reservoir cap.
- If you notice the presence of an air bubble at the level of the glass membrane or diaphragm, remove it by shaking the electrode.
- Connect the reference electrode and measurement electrode to the instrument.
- Open the opening tab on the reference electrode.
- Calibrate the electrodes by immersing them simultaneously in buffer solutions. Rinse them with demineralized water after each measurement.
- Immerse the tips of the electrodes in the solution to be measured, to the same height, while keeping them vertical and making sure that the membrane of the measurement electrode and the diaphragm of the reference electrode are immersed, but without exceeding the level of the internal electrolyte.
- Wait for the measurement to stabilize.
- Rinse the electrode with demineralized water after each measurement.
- Close the tab on the reference electrode again when you have finished measuring.

STATE AT DELIVERY

The separate pH measurement electrode is delivered in a cardboard box with:

- a shim,
- a storage reservoir cap,
- a bilingual user's manual,
- a multilingual accompanying datasheet.

For accessories and replacement parts, please visit our website:

www.chauvin-arnoux.com

CLEANING

The presence of white crystalline deposits on the electrode is a normal phenomenon due to evaporation of the storage solution in the reservoir cap. To remove this deposit, simply rinse the electrode in demineralized water.

Do not use aggressive and/or abrasive products (hydrofluoric acid, acetone, etc.). The glass membrane must not be scratched during cleaning.

Wear of the membrane is the most frequent cause of measurement problems.

If you encounter measurement problems, check that the membrane is clean. If there is a deposit, remove it. To do so, first rinse the electrode with water. If impurities are still present, use a solvent (such as ethanol) or a weak acid solution (diluted acetic acid or hydrochloric acid). It may sometimes be necessary to use a soft cloth but take care not to scratch the glass membrane.

MAINTENANCE

Keep your electrode in a vertical position in a reservoir cap containing an electrolytic solution (KCl 1 mol/L).

Do not store an electrode in dry conditions or in demineralized water. The membrane may be affected and the life span of the electrode may be reduced.

Electrode reference	Storage solution
BV41A/BV41H	KCl 1 mol/L storage solution
BV42A/BV42H	
XV41	
XV42	

If your electrode has been stored in dry conditions, reactivate the membrane. Immerse the electrode in a KCl 1mol/L solution or in a pH 4 buffer solution for a few hours.

Regularly check that your electrode is in good condition and clean (no dust, no scratches, etc.).

LIFE SPAN AND WARRANTY

pH measurement electrodes are wear parts. The life span of your electrode depends on the conditions in which it is used and the maintenance which you carry out on it.

The warranty does not apply if the glass is broken.

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